

TP500 T/PUMP® and Mul•T•Pads®

Animal Warming System

Help prevent body heat loss
and hypothermia during
surgery and recovery.

Surround your patient with
continuous warmth by
connecting 2 Mul•T•Pads
to one T/Pump.



 **GAYMAR®**

T/Pump®

Animals receive the benefit of warming with the portable, easy to use T/Pump System.

- Key operated controller helps prevent tampering
- Two over temperature safety devices are activated if needed
- If tipped-over, heater automatically shuts off and resets when pump is properly oriented
- Immediately lights up in event of thermal malfunction



Leakproof Connections

- Patented Klik-Tite® connectors snap into place quickly and easily



Secure Hose Attachments

- Permanently attached hoses assure proper, secure connections

T/Pump® Specifications

Dimensions: 8.12" W x 6.28" H x 5.62" D
(21 x 16 x 14cm)

Weight: 5lbs. 2oz. (2.3kg)

Hose Length: 10' (3.0m)

Power: 120 VAC, 60Hz, 200 Watts, 3 amp circuit breaker, 9' (2.7m) power cord (18/3 SJT3)

Fluid System: Distilled water, 51oz. reservoir (1500ml)
9-14 gph flow rate (34-56 lph flow rate)
Wet or dry storage
Water level indicator window

Control: ±2°F (1°C) accuracy

System: 85°F - 107°F (29.4°C - 41.7°C) heating range

Ratings: UL and ULC listed

T/Pump® Ordering Information

| MODEL NO. | DESCRIPTION |
|-----------|------------------------------|
| TP500 | with Klik-Tite connectors |
| TP500C | with Colder style connectors |
| TPP522 | TP500 with TP22C pad |
| TP472 | 230 VAC, 50Hz, "CE" marked |



Gaymar Industries, Inc.
10 Centre Drive • Orchard Park, NY 14127
Toll Free 1.800.828.7341
Fax 1.716.662.8624
www.gaymar.com

Mul•T•Pad®

As part of the T/Pump System, Mul•T•Pads® deliver a consistent flow of continuous warmth.

- The only pads available that can connect to another pad to surround your patient with effective warming (without special accessories when using Klik-Tite connectors)
- Pads are easily wiped down for multiple applications
- Unique flow design provides optimum warming even when folded



Mul•T•Pad® Ordering Information

| MODEL NO. | DESCRIPTION |
|--|---|
| Mul•T•Pads with Klik-Tite® connectors | |
| TP22B | 15" x 22" (38 x 56cm), Pliable All-polymer |
| TP22C | 15" x 22" (38 x 56cm), Pliable polymer on one side/ nonwoven fabric on the other |
| TP3E | 3" x 23" (8 x 58cm), Pliable polymer on one side/ nonwoven fabric on the other |
| TP12E | 13" x 18" (33 x 46cm), Pliable polymer on one side/ nonwoven fabric on the other |
| TP22E | 15" x 22" (38 x 56cm), Pliable polymer on one side/ nonwoven fabric on the other |
| TP26E | 18" x 26" (46 x 66cm), Pliable polymer on one side/ nonwoven fabric on the other |
| TP22G | 15" x 22" (38 x 56cm), Heavy duty pliable polymer |
| HP7020 | 25" x 33" (64 x 84cm), Heavy duty polymer |

Mul•T•Pads with Colder style connectors

| | |
|--------|---|
| TP612E | 13" x 18" (33 x 46cm), Nonwoven fabric both sides |
| TP622E | 15" x 22" (38 x 56cm), Nonwoven fabric both sides |
| TP626E | 18" x 26" (46 x 66cm), Nonwoven fabric both sides |

Accessories

| | |
|------|-----------------------------|
| TPS1 | Temperature therapy stand |
| HBX8 | 10' (3.05m) Hose attachment |

Ask about Gaymar's Convective Air
Animal Warming System!

Thermacare®



For availability in Canada, call:
Waterloo Bedding Co. Ltd.
825 Trillium Drive • Kitchener, ON N2R 1J9
Toll Free: 1.800.203.4293 • 519.742.4447
Fax: 519.742.3699