



**Meena Medical
Equipment, Inc.**

1-888-225-2502
www.meenamedical.com

CRITICARE, 5040 PONI MANUFACTURER SPECIFICATIONS

MODEL	504o PONI
ALERTS (CONFIGURATION)	
Audiovisual (ALERTS)	Hi/lo SpO2 and pulse rate, sensor off, sensor disconnect
Visual (ALERTS)	Audio off, low battery, hi/lo SpO2 and pulse rate, sensor off, sensor disconnect
Kind of battery (BATTERY)	Rechargeable, sealed lead-acid
Life, hours (ALERTS)	3-4 on full charge
DIMENSIONS (HXWXD)CM, (IN) (ANALYZER)	8x 13 x 18.5 (3.5x 5.5 x 7.5)
WEIGHT, kg (lb) (DISPLAY)	1.5(2.3) able;
OTHER SPECIFICATIONS (AIR SUPPLY)	12-trend memory;
DISPLAYS (STAINING POSTS)	SpO2, pulse rate and other parameters. LED, SpO2 and pulse rate; pulse amplitude bar; indicators for AC power, low battery, alarms off, battery use; system status; hi/lo arrows
SpO2 RANGE, % (KIND OF BATTERY)	0
Accuracy, % (POWER NEEDED)	70-99±2. 40-69±3
RESPONSE TIME, sec	Beat to beat
ALERT OVERRIDE (PROBES)	Yes
Reactivation (PROBES)	Automatic after 2 min or manual
VOLUME REGULATION (POWER NEEDED)	Yes
PROBES	Adult finger, Shell sensors, multisite reusable, disposable
Patients (PROBES)	Adult, pediatric, infant, neonate
Cable length, m (ft) (WORKING RADIUS,)	0.9(3),3.05(10), 4.5(15)
LINE POWER, VAC (TEMPERATURE)	110/120/220/240
EXTERNAL OUTPUT (PROBES)	Analog, digital RS232
PRINTER/RECORDER (CONFIGURATION)	No
SELF-TEST MODE (KIND OF BATTERY)	Yes
SETTLING TIME (PROBES)	4-6 sec from no pulse to numbers